



## TECAN Spark Plate Reader Configuration for Valita®Titer

Parameter	Instrument Settings	
Product	Valita®Titer (#VAL003)	Valita®Titer Plus (#VAL004)
<b>Plate type</b>	Corning/Costar 96-well Half Area Black Flat Bottom, #3694	
<b>Mode</b>	Fluorescence Polarization	Fluorescence Polarization
<b>Excitation</b>	Filter	Monochromator [optimum]; filter
<b>Excitation Wavelength</b>	485nm	485nm
<b>Excitation bandwidth</b>	20nm	20nm
<b>Emission</b>	Filter	Monochromator [optimum]; filter
<b>Emission Wavelength</b>	535nm	535nm
<b>Emission bandwidth</b>	20nm	5-20nm
<b>Flash Number</b>	>/= 100	200
<b>Integration</b>	40µs	40µs
<b>Gain</b>	Optimised using most fluorescent well [0mg/L]	Optimised using most fluorescent well [0mg/L]
<b>Blank</b>	Not defined	Not defined
<b>RFU [%]</b>	70	70
<b>Z position</b>	Optimised using most fluorescent well [0mg/L]	Optimised using most fluorescent well [0mg/L]
<b>Mirror</b>	Dichroic 510	Dichroic 510; 50% Mirror
<b>G-Factor</b>	<p>Manual: 1. If G-factor optimization is required, follow the instructions for Valita®TITER Plus (#VAL004).</p>	<p>Optimisation of G-factor required. Measure a blank well [120 µL of media in a Corning/Costar 96-well Half Area Black Flat Bottom, #3694] and a reference well [120 µL of media in Valita®TITER plate]. If a plate "Corning/Costar 96-well Half Area Black Flat Bottom, #3694" is not available Valita®TITER probe can be removed from Valita®TITER plate following the below procedure:</p> <ul style="list-style-type: none"> <li>• Add 200 µl of fresh media in a Valita®TITER plate well.</li> <li>• Mix it well with a pipette;</li> <li>• Incubate at room temperature for 30-minutes;</li> <li>• Mix it well with a pipette;</li> <li>• Wash with PBS (x3).</li> </ul> <p>Avoid shaking and/or 37°C incubation to avoid damaging the rest of the plate.</p>