



## Molecular Devices Paradigm Microplate Reader Configuration for Valita®Titer

Parameter	Instrument Settings	
Product	<b>Valita®Titer (#VAL003 / VAL013)</b>	<b>Valita®Titer Plus (#VAL004 / VAL014)</b>
Plate type	96-well: 96 Well Corning Half Area opaque 384 well: 384 Well Corning opaque	
Measurement Mode	Fluorescence Polarization	
Excitation	485 nm	
Emission	535 nm	
Integration time	400 ms	
G-Factor	1.0	
Read height (mm)	4.1 mm	4.8 mm
Blank correction	None; 0mg/L standard should be labelled as standard and not blank	